

IHE Integration Statement

<i>Vendor</i>	<i>Product Name</i>	<i>Version</i>
Fujifilm Medical Systems	FN-PS 551 Print Server	A07
<i>This product implements all of the transactions specified in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options specified below:</i>		
<i>Integration Profiles Implemented</i>	<i>Actors Implemented</i>	<i>Options Implemented</i>
Consistent Presentation of Images	Print Server	
<i>Links to Standards Conformance Statements of the Implementation</i>		
<i>HL7</i>		
<i>DICOM</i>	http://www.fujimed.com/sub/FN-PS551_DCS.pdf	
<i>Link to further information on this product:</i>	http://www.fujimed.com/medical/ldimagers.html	
<i>Date of Statement:</i>	November 1, 2002	