

## IHE Integration Statement

Vendor	Product Name	Version	Date
FUJIFILM Medical Systems U.S.A., Inc.	Synapse	4.x	11/18/2011
The product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Image Manager Image Display Performed Procedure Step Manager		
Patient Information Reconciliation	Image Manager Performed Procedure Step Manager		
Access to Radiology Information	Image Manager Report Reader		
Simple Image and Numeric Report	Report Reader		
Consistent Time	Time Client		
Teaching File and Clinical Trial Export	Export Selector		
Consistent Presentation of Images	Image Manager		
Key Image Note	Image Manager		
Mammography Image	Image Manager		
Portable Data for Imaging	Portable Media Creator		
<b>Internet address for vendor's IHE information:</b> <a href="http://www.fujifilmusa.com/products/medical/technical-documents/product-bulletins/">http://www.fujifilmusa.com/products/medical/technical-documents/product-bulletins/</a>			
<b>Links to Standards Conformance Statements of the Implementation</b>			
DICOM	<a href="http://www.fujifilmusa.com/products/medical/technical-documents/product-bulletins/">http://www.fujifilmusa.com/products/medical/technical-documents/product-bulletins/</a>		
HL7	<a href="http://www.fujifilmusa.com/products/medical/technical-documents/product-bulletins/">http://www.fujifilmusa.com/products/medical/technical-documents/product-bulletins/</a>		
Links to general information on IHE			
In North America: <a href="http://www.ihe.net">www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.ihe-j.org">www.ihe-j.org</a>	