

| <b>IHE Integration Statement</b>  |   |  |             |
|---|---|--|-------------|
| <b>Vendor</b>   | <b>Product Name</b>   | <b>Version</b>   | <b>Date</b> |
| FUJIFILM Corporation  | FDX Console / Console Advance<br>(DR-ID 300CL)  | 5.0  | 5/26/2011   |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:              |   |  |             |
| <b>Integration Profiles Implemented</b>   | <b>Actors Implemented</b>   | <b>Options Implemented</b>   |             |
| Scheduled Workflow (SWF)  | Acquisition Modality  | Broad Worklist Query<br>Modality Group Case                                      |             |
| Patient Information Reconciliation (PIR)  |   |  |             |
| <b>Internet address for vendor's IHE information:</b> <a href="http://www.fujimed.com/sub/product_bulletins.html">http://www.fujimed.com/sub/product_bulletins.html</a> |   |  |             |
| <b>Links to Standards Conformance Statements for the Implementation</b>   |   |  |             |
| <b>HL7</b>  |   |  |             |
| <b>DICOM</b>  | <a href="http://www.fujimed.com/sub/product_bulletins.html">http://www.fujimed.com/sub/product_bulletins.html</a> |  |             |
| <b>Links to general information on IHE</b>  |   |  |             |
| In North America: <a href="http://www.ihe.net">www.ihe.net</a>  | In Europe: <a href="http://www.ihe-europe.net">www.ihe-europe.net</a>   | In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a> |             |