

## IHE Integration Statement

Vendor	Product Name	Version	Date
FUJIFILM Medical Systems USA, Inc.	Synapse	3.2.1 SR-334	24 March 2009

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow	Image Manager Image Display Performed Procedure Step Manager	
Patient Information Reconciliation	Image Manager Performed Procedure Step Manager	
Access to Radiology Information	Image Manager Report Reader	
Simple Image and Numeric Report	Report Reader	
Consistent Time	Time Client	
Teaching File and Clinical Trial Export	Export Selector	
Consistent Presentation of Images	Image Manager	
Key Image Note	Image Manager	
Mammography Image	Image Manager	
Portable Data for Imaging	Portable Media Creator	

**Internet address for vendor's IHE information:** [http://www.fujimed.com/sub/product\\_bulletins.html](http://www.fujimed.com/sub/product_bulletins.html)

### Links to Standards Conformance Statements of the Implementation

<b>DICOM</b>	<a href="http://www.fujimed.com/sub/product_bulletins.html">http://www.fujimed.com/sub/product_bulletins.html</a>
--------------	---

### Links to general information on IHE

In North America: <a href="http://www.rsna.org/IHE">www.rsna.org/IHE</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>
--	---	--