

Fujifilm - IHE Integration Statement		Date
		07 March 2008
Vendor	Product Name	Version
Fujifilm Medical Systems	Synapse	Synapse 3.2.1
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow	Image Manager Image Display Performed Procedure Step Manager	
Patient Information Reconciliation	Image Manager Performed Procedure Step Manager	
Access to Radiology Information	Image Manager Report Reader	
Simple Image and Numeric Report	Report Reader	
Consistent Time	Time Client	
Teaching File and Clinical Trial Export	Export Selector	
Mammography Image	Image Manager	
PDI	Portable Media Creator	
Links to Standards Conformance Statements of the Implementation		
HL7	http://www.fujimed.com/	
DICOM	http://www.fujimed.com/	
Links to General Information on IHE		
In North America www.rsna.org/IHE		
In Europe: www.ihe-europe.org		
In Japan: www.jira-net.or.jp/ihe		